

Canine CRP



Application Note for the Gentian Canine CRP Immunoassay on the Alinity¹

For *in vitro* diagnostic use by laboratory professionals.

This document describes the instrument specific settings and performance of the product on the instrument above. For assay information, please refer to the IFU available on www.gentian.com.

Assay kit components

Products available	
Gentian Canine CRP Reagent Kit <ul style="list-style-type: none">R1 Assay Buffer (45 mL)R2 Immunoparticles (10.5 mL)	REF 1501
Gentian Canine CRP Calibrator Kit (6 levels x 0.5 mL)	REF 1551
Gentian Canine CRP Control Kit (2 levels x 0.5 mL)	REF 1519

All products are ready for use.

Reagent stability

The in-use stability of the Gentian Canine CRP Reagent Kit was found to be at least 4 weeks in an on board study based on the CLSI guideline EP25 [1].

Calibration stability

The calibration curve stability of the Gentian Canine CRP Calibrator Kit was found to be at least 2 weeks in a study based on the CLSI guideline EP25 [1].

Performance characteristics

All results refer to validation of the Gentian Canine CRP Immunoassay on one instrument site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Canine CRP Immunoassay was found to be 9-292 mg/L. The exact measuring range is specific to the calibrator lot, please refer to the analytical value sheet available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian Canine CRP Immunoassay was tested in a study based on the CLSI guideline EP17 [2]. The limit of quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian Canine CRP Immunoassay was found to be 8 mg/L.

Linearity

The linearity range of the Gentian Canine CRP Immunoassay was found to be 9-292 mg/L in a linearity study based on the CLSI guideline EP06 [3].

Security zone

No antigen excess effect in samples below 913 mg/L was observed for the Gentian Canine CRP Immunoassay in a study based on the CLSI guideline EP34 [4]. Samples with a CRP concentration above the highest calibrator and up to 913 mg/L return a value above the highest calibrator and are flagged for rerun.

Precision

Precision of the Gentian Canine CRP Immunoassay was tested in a 3-day precision study based on the CLSI guideline EP05 [5]. 3 serum pools and 2 controls were measured 5 times with 5 replicates (n=25).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
S1	13	4	4	5
S2	37	1	3	3
S3	184	0	2	2
CL	30	1	0	1
CH	100	1	2	2

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [6]. The Gentian Canine CRP Immunoassay had a recovery of 86-102 %.

Analytical specificity and limitations

Interference was tested in a study based on the CLSI guideline EP07 [7]. As the antibodies in the Gentian Canine CRP Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [8]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	5 g/L
Intralipid	10 g/L
Bilirubin	0.6 g/L

Instrument variation

Results obtained with the Gentian Canine CRP Immunoassay were compared using Passing-Bablok regression with results from the Cobas c501 instrument (Roche) in a study based on the CLSI guideline EP09 [9].

n	Range of samples [mg/L]	Term	Coefficient	95% CI
40	9.7-294.1	Intercept	-0.31	[-1.01, 0.45]
		Slope	1.04	[1.03, 1.06]
		R ²	1.00	



Bjornasveien 5
N-1596 Moss
Norway
TEL: +47 99 33 99 05
www.gentian.com

References

1. CLSI. Evaluation of Stability of *In Vitro* Diagnostic Reagents; Approved Guideline. CLSI document EP25-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
2. CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012
3. CLSI. Evaluation of Linearity of Quantitative Measurement Procedures. 2nd ed. CLSI guideline EP06. Clinical and Laboratory Standards Institute; 2020
4. CLSI. Establishing and verifying an extended measuring interval through specimen dilution and spiking. 1st ed. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
5. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014
6. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258
7. CLSI. Interference Testing in Clinical Chemistry. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
8. Larsson A, et al. Poultry Science 1993;72:1807-12
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

Modification from the previous version

- First version.

Date of issue

2025-03-19

Instrument Settings for the Gentian Canine CRP Immunoassay on the Alinity¹

GENERAL	CALIBRATION	RESULTS	RETEST RULES	SMART WASH														
Assay Name: ** Assay Number: *																		
Assay Type: Photometric Assay Version: 1																		
Assay Status: Primary Result Units: mg/L																		
Date/Time: * Operator: *																		
Run Controls for Onboard Reagents by: **																		
Assay Availability: Enable																		
Reaction Definition																		
Reaction Mode: End Up Absorbance Range: - Sample Blank Types: Self Blank																		
Primary Wavelength: 640 Last Read: 31 Blank Read Time: 18																		
Secondary Read Times Main: 30 - 31 Flex: - Color Correction: -																		
Reagent																		
Reagent: ** Diluent: None Diluent Dispense Mode: **																		
Validity Checks																		
Reaction Check: None A Read Time: - B Read Time: -																		
Calibration Limits:																		
Minimum: Maximum:																		
Sample																		
<table border="1"> <thead> <tr> <th>Dilution Name</th> <th>Sample</th> <th>Diluted sample</th> <th>Diluent</th> <th>Water</th> <th>Dilution Factor</th> <th>Default Dilution</th> </tr> </thead> <tbody> <tr> <td>BSA</td> <td>2,0</td> <td></td> <td></td> <td></td> <td>1,0</td> <td>-</td> </tr> </tbody> </table>					Dilution Name	Sample	Diluted sample	Diluent	Water	Dilution Factor	Default Dilution	BSA	2,0				1,0	-
Dilution Name	Sample	Diluted sample	Diluent	Water	Dilution Factor	Default Dilution												
BSA	2,0				1,0	-												
Reagent Volume: 250 Water Volume: - Dispense Mode: Type 1																		
R1 R2 Type 1																		

GENERAL	CALIBRATION	RESULTS	RETEST RULES	SMART WASH
New Reagent				
Reagent Name: ** Reagent Type: R1 and R2 Low Alert: ** Number of Tests: 100 Onboard Stability: 1200 Run Calibration for Reagent by: **				

* Default by instrument

** User defined

*** Lot specific. See analytical value sheet available on www.gentian.com

GENERAL	CALIBRATION	RESULTS	RETEST RULES	SMART WASH																																								
Assay Name: ** Assay Number: *	Assay Type: Photometric Assay Version: *	Assay Status: Primary Result units: mg/L	Date/Time: * Operator: *																																									
Calibration																																												
Calibration Method: Spline Factor:		Full Calibration Interval: 999 (Hours) Default Ordering Type:																																										
Use Cal Factor From: Calibration Expiration Warning: ** (Hours)		Adjust Type: None Adjust Interval: (Hours) Adjust Level:																																										
Calibrators		Validity Checks																																										
Calibrator Set: ** Replicates: 2		Blank Absorbance Range: - Span: Blank - Span absorbance range: - Expected Cal Factor: Expected Cal Factor Tolerance %: Maximum Curve Fit:																																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 20%;">Calibrator Level:</th> <th style="width: 10%;">Sample</th> <th style="width: 10%;">Diluted Sample</th> <th style="width: 10%;">Diluent</th> <th style="width: 10%;">Water</th> </tr> </thead> <tbody> <tr> <td>Blank:</td> <td>***</td> <td>2,0</td> <td></td> <td></td> </tr> <tr> <td>Cal 1:</td> <td>***</td> <td>2,0</td> <td></td> <td></td> </tr> <tr> <td>Cal 2:</td> <td>***</td> <td>2,0</td> <td></td> <td></td> </tr> <tr> <td>Cal 3:</td> <td>***</td> <td>2,0</td> <td></td> <td></td> </tr> <tr> <td>Cal 4:</td> <td>***</td> <td>2,0</td> <td></td> <td></td> </tr> <tr> <td>Cal 5:</td> <td>***</td> <td>2,0</td> <td></td> <td></td> </tr> <tr> <td>Cal 6:</td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>					Calibrator Level:	Sample	Diluted Sample	Diluent	Water	Blank:	***	2,0			Cal 1:	***	2,0			Cal 2:	***	2,0			Cal 3:	***	2,0			Cal 4:	***	2,0			Cal 5:	***	2,0			Cal 6:				
Calibrator Level:	Sample	Diluted Sample	Diluent	Water																																								
Blank:	***	2,0																																										
Cal 1:	***	2,0																																										
Cal 2:	***	2,0																																										
Cal 3:	***	2,0																																										
Cal 4:	***	2,0																																										
Cal 5:	***	2,0																																										
Cal 6:																																												

GENERAL	CALIBRATION	RESULTS	RETEST RULES	SMART WASH
Assay Name: ** Assay Number: *	Assay Type: Photometric Assay Version: *	Assay Status: Primary Result units: mg/L	Date/Time: * Operator: *	
Results Unit				
Results Units: mg/L Correlation factor: 1,0000		Decimal Places: 1 Intercept: 0,0000		Results Unit UCUM:
Results				
Gender	Age Unit:	Low Linearity: Age Range	Normal Range	High Linearity: Extreme Range
Interpretation				
Name		Range	Review Required	
<NotDefined>		< 0,00	-	
<NotDefined>		>= 0,00	-	
<NotDefined>		>= 0,00	-	
<NotDefined>		>= 0,00	-	
<NotDefined>		>= 0,00	-	

* Default by instrument

** User defined

*** Lot specific. See analytical value sheet available on www.gentian.com

Disclaimer: The specific settings above is what used to validate the application on the specific instrument. For any instrument specific settings, please refer to the instrument manual. Please be aware that illustrations or settings might be affected in case of an instrument software update.