

# Application Note for the Gentian Calprotectin Immunoassay on the AU680<sup>1</sup>

For *in vitro* diagnostic use by laboratory professionals.

This document describes the instrument specific settings and performance of the product on the instrument above. For assay information, please refer to the IFU available on [www.gentian.com](http://www.gentian.com).

## Assay kit components

Products available	
Gentian GCAL® Calprotectin Reagent Kit	REF 1201
<ul style="list-style-type: none"> <li>R1 Assay Buffer (54 mL)</li> <li>R2 Immunoparticles (9 mL)</li> </ul>	
Gentian GCAL® Calprotectin Reagent Kit S	REF 1202
<ul style="list-style-type: none"> <li>R1 Assay Buffer (30 mL)</li> <li>R2 Immunoparticles (5 mL)</li> </ul>	
Gentian GCAL® Calprotectin Calibrator Kit (6 levels x 1 mL)	REF 1251
Gentian GCAL® Calprotectin Control Kit (2 levels x 1 mL)	REF 1219

All products are ready for use.

## Reagent stability

The in-use stability of the Gentian GCAL® Calprotectin Reagent Kit was found to be at least 4 weeks in an on-board study based on the CLSI guideline EP25 [1]. If the instrument remains unused for more than a week, please ensure the reagents are gently inverted every 7 days.

## Calibration stability

The calibration curve stability of the Gentian GCAL® Calprotectin Calibrator Kit was found to be at least 1 week in a study based on the CLSI guideline EP25 [1].

## Performance characteristics

All results refer to validation of the Gentian GCAL® Calprotectin Immunoassay on one instrument site with one lot of reagents, unless otherwise stated.

## Measuring range

The measuring range of the Gentian GCAL® Calprotectin Immunoassay was found to be 0.48-19.16 mg/L. The exact measuring range is specific to the calibrator lot, please refer to the analytical value sheet available on [www.gentian.com](http://www.gentian.com).

## Analytical sensitivity

The analytical sensitivity of the Gentian GCAL® Calprotectin Immunoassay was tested in a study based on the CLSI guideline EP17 [2]. The limit of quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian GCAL® Calprotectin Immunoassay was found to be 0.35 mg/L for lithium heparin plasma and serum.

## Linearity

The linearity range of the Gentian GCAL® Calprotectin Immunoassay was found to be 0.44-22.02 mg/L for lithium heparin plasma and 0.41-20.98 mg/L for serum in a linearity study based on the CLSI guideline EP06 [3].

## Security zone

No antigen excess effect in samples below 90 mg/L was observed for the Gentian GCAL® Calprotectin Immunoassay in a study based on the CLSI guideline EP34 [4]. Samples with a calprotectin concentration above the highest calibrator and up to 90 mg/L return a value above the highest calibrator and are flagged for rerun.

## Precision

Precision of the Gentian GCAL® Calprotectin Immunoassay was tested in a 3-day precision study based on the CLSI guideline EP05 [5]. 3 lithium heparin plasma (P1, P2 and P3) and 3 serum (S1, S2 and S3) pools and 2 controls (CL, CH) were measured 5 times with 5 replicates (n=25).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.90	3.50	1.27	3.72
P2	6.01	0.67	0.97	1.18
P3	13.64	0.71	1.87	2.00
S1	1.02	3.57	0.00	3.57
S2	5.84	0.61	0.73	0.95
S3	14.27	0.86	1.53	1.76
CL	1.01	3.23	3.00	4.41
CH	9.98	0.81	0.57	0.99

## Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [6]. The Gentian GCAL® Calprotectin Immunoassay had a recovery of 95-97 %.

## Analytical specificity and limitations

Interference was tested in a study based on the CLSI guideline EP07 [7]. As the antibodies in the Gentian GCAL® Calprotectin Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [8]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	2.5 g/L
Intralipid	10 g/L
Bilirubin	0.6 g/L

## Instrument variation

Results obtained with the Gentian GCAL® Calprotectin Immunoassay were compared using Passing-Bablok regression with results from the Cobas c501 instrument (Roche) in a study based on the CLSI guideline EP09 [9].

n	Range of samples [mg/L]	Term	Coefficient	95% CI
50	0.50-21.07	Intercept	-0.12	[-0.19, -0.05]
		Slope	1.13	[1.12, 1.15]
		R <sup>2</sup>	1.00	

Bjornasveien 5  
N-1596 Moss  
Norway  
TEL: +47 99 33 99 05  
[www.gentian.com](http://www.gentian.com)

## References

1. CLSI. Evaluation of Stability of *In Vitro* Diagnostic Reagents; Approved Guideline. CLSI document EP25-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
2. CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012
3. CLSI. Evaluation of Linearity of Quantitative Measurement Procedures. 2<sup>nd</sup> ed. CLSI guideline EP06. Clinical and Laboratory Standards Institute; 2020
4. CLSI. Establishing and verifying an extended measuring interval through specimen dilution and spiking. 1<sup>st</sup> ed. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
5. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014
6. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258
7. CLSI. Interference Testing in Clinical Chemistry. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
8. Larsson A, et al. Poultry Science 1993;72:1807-12
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

## Modification from the previous version

- Harmonised analytical measuring range across Beckman Coulter instruments

## Date of issue

2024-11-26

## Instrument Settings for the Gentian GCAL® Calprotectin Immunoassay on the AU680<sup>1</sup> (serum/plasma)

Reagent ID: 254

Specific Test Parameters											
General		LIH	ISE	Range							
Test Name:		CAL1G	<	>	Type:		Serum***	Operation:		Yes	
Sample Volume		4.0	μL	Dilution		0	μL	OD Limit			
Pre-Dilution Rate		1						Min. OD		Max. OD	
Reagents Volume:		R1(R1-1)	180	μL	Dilution		0	μL	Reagent OD limit:		
								First Low		-2.0000 High 2.5000	
								Last Low		-2.0000 High 2.5000	
R2 Volume		30	μL	Dilution		0	μL	Dynamic Range Low		0.48 High 19.16	
Common Reagent		Type	None	Name			Correlation Factor A		1 B 0		
Wavelength:		Pri.	660	nm	Sec.		None	nm	Factor for Maker A 1 B 0		
Method:		END									
Reaction slope:		+				Onboard Stability		28	Days	0	Hour
Measuring Point 1:		First	9	Last		15	LIH Influence Check		No		
Measuring Point 2:		First		Last			Lipemia				
Linearity:						Icterus					
No Lag Time:						Hemolysis					

Specific Test Parameters											
General		ISE		Range							
Test Name:		CAL1G ▾		<div style="display: inline-block; width: 30px; text-align: center;">&lt;</div> <div style="display: inline-block; width: 30px; text-align: center;">&gt;</div>		Type:		Serum*** ▾			
Value/Flag:		* ▾		Level L:		*		Level H:		*	
Specific Ranges:											
		From		To							
	Sex	Year	Month	Year	Month	Low	High				
<input type="checkbox"/>	1.	* ▾	* ▾	* ▾	* ▾	* ▾	* ▾				
<input type="checkbox"/>	2.	* ▾	* ▾	* ▾	* ▾	* ▾	* ▾				
<input type="checkbox"/>	3.	* ▾	* ▾	* ▾	* ▾	* ▾	* ▾				
<input type="checkbox"/>	4.	* ▾	* ▾	* ▾	* ▾	* ▾	* ▾				
<input type="checkbox"/>	5.	* ▾	* ▾	* ▾	* ▾	* ▾	* ▾				
<input type="checkbox"/>	6.	* ▾	* ▾	* ▾	* ▾	* ▾	* ▾				
7. No demographics						* ▾	* ▾				
8. Not within expected values						* ▾	* ▾				
Unit		mg/L		Decimal Places		*					

Panic Value

Low

\*

High

\*

Calibration Specific											
General		ISE									
Test Name:		CAL1G ▾		< >		Type		Serum*** ▾		<input type="checkbox"/> Use Serum Cal.	
Calibration Type:		6AB ▾		Formula:		Spline ▾		Counts:		2	
<Calibrator Parameters>											
	Calibrator	OD	Conc	Factor Range				Slope Check		None ▾	
				Low	High						
Point 1:	1 ▾		**	-2.0000	2.5000			Allowable Range Check			
Point 2:	2 ▾		**	-2.0000	2.5000			<input type="checkbox"/> Reagent Blank			
Point 3:	3 ▾		**	-2.0000	2.5000			<input type="checkbox"/> Calibration			
Point 4:	4 ▾		**	-2.0000	2.5000			Advanced Calibration			
Point 5:	5 ▾		**	-2.0000	2.5000			Operation		* ▾	
Point 6:	6 ▾		**	-2.0000	2.5000			Interval (RB/ACAL)		* ▾	
Point 7:	▾										
Point 8:	▾										
Point 9:	▾										
Point 10:	▾										
<Point Cal. For Master Curve>											
Calibrator		OD	Conc	OD Range		Use Master Curve		Stability		<input type="checkbox"/> Lot Calibration	
				Low	High						
Point 1:								Reagent Blanks		7 Day 0 Hour	
Point 2:								Calibration		7 Day 0 Hour	
MB Type Factor:			1-Point Calibration Point		None ▾		<input type="checkbox"/> With CONC-0				

*Disclaimer: The specific settings above is what used to validate the application on the specific instrument. For any instrument specific settings, please refer to the instrument manual. Please be aware that illustrations or settings might be affected in case of an instrument software update.*

\* User defined

\*\* Lot specific. See analytical value sheet available on [www.gentian.com](http://www.gentian.com)

\*\*\*Valid for both serum and lithium heparin plasma