

Application Note for the Gentian Calprotectin Immunoassay on the DxC 500 AU¹

For *in vitro* diagnostic use by laboratory professionals.

This document describes the instrument specific settings and performance of the product on the instrument above. For assay information, please refer to the IFU available on www.gentian.com.

Assay kit components

Products available	
Gentian GCAL [®] Calprotectin Reagent Kit	REF 1201
<ul style="list-style-type: none"> R1 Assay Buffer (54 mL) R2 Immunoparticles (9 mL) 	
Gentian GCAL [®] Calprotectin Reagent Kit S	REF 1202
<ul style="list-style-type: none"> R1 Assay Buffer (30 mL) R2 Immunoparticles (5 mL) 	
Gentian GCAL [®] Calprotectin Calibrator Kit (6 levels x 1 mL)	REF 1251
Gentian GCAL [®] Calprotectin Control Kit (2 levels x 1 mL)	REF 1219

All products are ready for use.

Reagent stability

The in-use stability of the Gentian GCAL[®] Calprotectin Reagent Kit was found to be at least 4 weeks in an on board study based on the CLSI guideline EP25 [1]. If the instrument remains unused for more than a week, please ensure the reagents are gently inverted every 7 days.

Calibration stability

The calibration curve stability of the Gentian GCAL[®] Calprotectin Calibrator Kit was found to be at least 1 week in a study based on the CLSI guideline EP25 [1].

Performance characteristics

All results refer to validation of the Gentian GCAL[®] Calprotectin Immunoassay on one instrument site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian GCAL[®] Calprotectin Immunoassay was found to be 0.48-19.16 mg/L. The exact measuring range is specific to the calibrator lot, please refer to the analytical value sheet available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian GCAL[®] Calprotectin Immunoassay was tested in a study based on the CLSI guideline EP17 [2]. The limit of quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian GCAL[®] Calprotectin Immunoassay was found to be 0.45 mg/L for lithium heparin plasma and 0.43 mg/L for serum.

Linearity

The linearity range of the Gentian GCAL[®] Calprotectin Immunoassay was found to be 0.41-21.39 mg/L in lithium heparin plasma and 0.38-20.25 mg/L in serum in a linearity study based on the CLSI guideline EP06 [3].

Security zone

No antigen excess effect in samples below 88 mg/L was observed for the Gentian GCAL[®] Calprotectin Immunoassay in a study based on the CLSI guideline EP34 [4]. Samples with a calprotectin concentration above the highest calibrator and up to 88 mg/L return a value above the highest calibrator and are flagged for rerun with automatic dilution.

Precision

Precision of the Gentian GCAL[®] Calprotectin Immunoassay was tested in a 3-day precision study based on the CLSI guideline EP05 [5]. 2 lithium heparin plasma (P1 and P1) and 2 serum (S1 and S2) pools and 2 controls (CL, CH) were measured 5 times with 5 replicates (n=25).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
PR-S1	0.90	4.74	1.10	4.86
PR-S2	8.21	0.82	1.00	1.29
PR-P1	5.18	1.15	0.98	1.51
PR-P2	13.62	1.23	0.54	1.35
CL	0.99	3.21	0.00	3.21
CH	9.98	0.74	0.22	0.77

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [6]. The Gentian GCAL[®] Calprotectin Immunoassay had a recovery of 101-114 %.

Analytical specificity and limitations

Interference was tested in a study based on the CLSI guideline EP07 [7]. As the antibodies in the Gentian GCAL[®] Calprotectin Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [8]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	2.5 g/L
Intralipid	10 g/L
Bilirubin	0.6 g/L

Instrument variation

Results obtained with the Gentian GCAL[®] Calprotectin Immunoassay were compared using Passing-Bablok regression with results from the Cobas c501 instrument (Roche) in a study based on the CLSI guideline EP09 [9].

n	Range of samples [mg/L]	Term	Coefficient	95% CI
50	0.42-20.07	Intercept	0.01	[-0.08, 0.07]
		Slope	1.10	[1.09, 1.13]
		R ²	1.00	

Bjornasveien 5
N-1596 Moss
Norway
TEL: +47 99 33 99 05
www.gentian.com

References

1. CLSI. Evaluation of Stability of *In Vitro* Diagnostic Reagents; Approved Guideline. CLSI document EP25-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
2. CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012
3. CLSI. Evaluation of Linearity of Quantitative Measurement Procedures. 2nd ed. CLSI guideline EP06. Clinical and Laboratory Standards Institute; 2020
4. CLSI. Establishing and verifying an extended measuring interval through specimen dilution and spiking. 1st ed. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
5. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014
6. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258
7. CLSI. Interference Testing in Clinical Chemistry. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
8. Larsson A, et al. Poultry Science 1993;72:1807-12
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

Modification from the previous version

- Harmonised analytical measuring range across Beckman Coulter instruments

Date of issue

2024-11-26

Instrument Settings for the Gentian GCAL[®] Calprotectin Immunoassay on the DxC 500 AU¹ (serum/plasma)

Reagent ID: 254

TEST CONFIGURATION & CHEMISTRY DETAILS									
Assay Name	Test	Rev	Discipline		Chemistry				
Test ID	GCAL*		Calculated Result		<input type="checkbox"/>				
LIS Code	GCAL*		Result Type		Quantitative ▼				
UNITS AND RANGE SETTINGS									
Use Settings from	Serum*** ▼	Units	mg/L ▼	Decimal Places	x.xx ▼	Serum***			
Test Kind	General ▼	Revision	01						
Reagent Name	GCAL*	Reagent ID	254						
ABB Name	CAL1G	Parameter Long Name	Calprotectin C97497 CAL1G serum						
GENERAL PARAMETERS									
SAMPLE VOLUME		Sample Volume	3.0 μL	Dilution	0 ▼ μL	REACTION OD LIMIT			
REAGENT VOLUME		Predilution Rate	1 ▼			Low -2.0000 High 3.0000			
		R1-1	150 μL	Dilution	0 μL	REACTION BLANK OD LIMIT			
		R2-1	22 μL	Dilution	0 μL	First: Low -2.0000 High 2.5000			
						Last: Low -2.0000 High 2.5000			
WAVELENGTH		Primary	660 nm	Secondary	None nm	ANALYTICAL MEASURING RANGE			
						Low 0.48 High 19.16			
METHOD		END ▼				MANUFACTURER FACTOR			
						A 1 B 0			
REACTION SLOPE		+				REAGENT ONBOARD STABILITY			
						28 Days 0 Hours			
MEASURING POINT		Point 1: First 9		Last 15		LIH INFLUENCE CHECK			
		Point 2: First				<input type="checkbox"/> Perform LIH check			
Linearity Limit						Lipemia ▼			
						Icterus ▼			
						Hemolysis ▼			
Lag Time Check						<input type="checkbox"/> Perform Lag Time Check			
CALIBRATION PARAMETERS									
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4
mg/L	2	None	0	None	0	None	0	None	0
CALIBRATOR SPECIFIC					CALIBRATION OD AND CONCENTRATION PARAMETERS				
Calibration Type 6AB					Counts 2 ▼				
Formula Spline					MB Factor				
Calibrator Name					Positive Cutoff				
Add *									
<input type="checkbox"/> SLOPE CHECK					Number of Levels 6				
Slope Check None									
STABILITY AND INTERVAL									
Reagent Blank Stability 7 Days 0 Hours					Interval * ▼				
Calibration Stability 7 Days 0 Hours					Interval * ▼				
					OD DELTA CHECK				
					<input type="checkbox"/> Reagent Blank				
					<input type="checkbox"/> Calibration				
PROZONE CHECK PARAMETERS									
<input type="checkbox"/> Logic Check 1			<input type="checkbox"/> Logic Check 2			<input type="checkbox"/> Logic Check 3			
Check Points			Check Points			Check Points			
Point 1 *			Point 1 *			Point 1 *			
Point 2 *			Interval *			Interval *			
Point 3 *						Value 1 *			
						Value 2 *			
Limit Points			Limit Points			Limit Points			
Limit 1 *			Limit 1 *			Limit 1 *			
Limit 2 *			Limit 2 *			Limit 2 *			
Check Pattern									
Pattern *									

Disclaimer: The specific settings shown below is what used to validate the application on the specific instrument. For any instrument specific settings, please refer to the instrument manual. Please be aware that illustrations or settings might be affected in case of an instrument software update.

* User defined

** Lot specific. See analytical value sheet available on www.gentian.com

***Valid for both serum and lithium heparin plasma