## Gentian Cystatin C Immunoassay na sistemima AU (AU5800, AU680, AU480, DxC 500 AU, DxC 700 AU) kompanije Beckman Coulter®

REF B08179

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2023-10-12

Za druge jezike posetite:

[www.gentian.com/products/ifu/cystatin-c/beckmancoulter](http://www.gentian.com/products/ifu/cystatin-c/beckmancoulter)

## Predstavnici

Odgovorna osoba u UK  
Emergo Consulting (UK) Limited  
c/o Cr360 – UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
Ujedinjeno Kraljevstvo

CH REP

MedEnvoy Switzerland  
Gothardstrasse 28  
6302 Zug  
Švajcarska

## Literatura

1. Abrahamson M et al: Biochem J 1990;268:287-94.
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.
4. Filler G et al: Clin Biochem 2005 ;38 :1-8.
5. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567.
6. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011.
7. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute ;2012.
8. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures*. 2<sup>nd</sup> ed. CLSI guideline EP06. Clinical Laboratory Standards Institut;2020.
9. CLSI. *Establishing and Verifying and Extended Measuring Interval Through Specimen Dilution and Spiking*; 1<sup>st</sup> Edition. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Institute; 2018.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258.
12. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
13. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8.
14. CLSI. *Interference Testing in Clinical Chemistry*. 3<sup>rd</sup> ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
15. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3<sup>rd</sup> ed. CLSI guideline EP09c. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
16. Qiu X et al.: Oncotarget. 2017;8(42):72985-72999.
17. CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.
18. Nitsch D, et al. Am J Kidney Dis. Jun 2011;57(6):863-72.

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Sistemski reagens: B08179 ID reagensa: 228

TEST CONFIGURATION & CHEMISTRY DETAILS												
Assay Name	Test	Rev	Discipline	Chemistry								
Test ID	CYS		Calculated Result	<input type="checkbox"/>								
LIS Code	CYS											
UNITS AND RANGE SETTINGS												
Use Settings from	Serum	Units	mg/L	Decimal Places	x.xx	Plasma						
Test Kind	General	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch								
Reagent Name	CYS	Reagent ID	228	<input type="checkbox"/> FSE Test								
ABB Name	CYS1G	Parameter Long Name	Cystatin C B08179 CYS1G CYSC Serum									
Region	<input checked="" type="checkbox"/> US <input checked="" type="checkbox"/> OUS <input checked="" type="checkbox"/> AP <input type="checkbox"/> JP <input checked="" type="checkbox"/> EU <input type="checkbox"/> Other											
GENERAL PARAMETERS												
SAMPLE VOLUME			Sample Volume	2.0	µL	Dilution	0	µL	REACTION OD LIMIT			
REAGENT VOLUME			Predilution Rate	1				Low		High		
WAVELENGTH			R1-1	150	µL	Dilution	0	µL	REACTION BLANK OD LIMIT			
METHOD			R2-1	30	µL	Dilution	10	µL	First: Low	-2.0000	High	2.0000
REACTION SLOPE									Last: Low	-2.0000	High	2.0000
MEASURING POINT									ANALYTICAL MEASURING RANGE			
LINEARITY LIMIT									Low	0.38	High	7.84
LAG TIME CHECK									MANUFACTURER FACTOR			
PERFORM LIH CHECK									A	1	B	0
PERFORM LAG TIME CHECK									REAGENT ONBOARD STABILITY			
PERFORM LAG TIME CHECK									60** Days			0** Hours
PERFORM LAG TIME CHECK									LIH INFLUENCE CHECK			
PERFORM LAG TIME CHECK									<input type="checkbox"/> Perform LIH check			
PERFORM LAG TIME CHECK									Lipemia			+ ▼
PERFORM LAG TIME CHECK									Icterus			+ ▼
PERFORM LAG TIME CHECK									Hemolysis			+ ▼
CALIBRATION PARAMETERS												
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4			
mg/L	2	None	0	None	0	None	0	None	0			
CALIBRATOR SPECIFIC					CALIBRATION OD AND CONCENTRATION PARAMETERS							
Calibration Type					6AB	Counts	2	<input type="checkbox"/> Use highest calibrator for Upper AMR				
Formula					Spline	MB Factor						
Calibrator Name					Add	CYS	Positive Cutoff					
<input checked="" type="checkbox"/> SLOPE CHECK					Number of Levels	6						
Slope Check												
STABILITY AND INTERVAL												
Reagent Blank Stability					28	Days	0	Hours	Interval	Lot		
Calibration Stability					28	Days	0	Hours	Interval	Lot		
OD DELTA CHECK												
<input type="checkbox"/> Reagent Blank							0.0000					
<input type="checkbox"/> Calibration							0.0000					
PROZONE CHECK PARAMETERS												
<input type="checkbox"/> Logic Check 1			<input type="checkbox"/> Logic Check 2			<input type="checkbox"/> Logic Check 3						
Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values			
Point 1	#	Value 1	#	Point 1	#	Value 1	#	Point 1	#			
Point 2	#	Value 2	#	Interval	#	Value 2	#	Interval	#			
Point 3	#	Value 3	#	Limit Points	Limit 1	#	Limit 1	#	Limit 1			
Limit Points	Limit 1	#	Limit 2	#	Limit 2	#	Limit 2	#	Limit 2			
Check Pattern	Pattern	#										

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Sistemski reagens: B08179 ID reagensa: 228

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> CYS1G <input type="button" value="Test No"/> <b>Type:</b> Serum <b>Operation:</b> Yes				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit <sup>1st</sup> Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.40"/> High <input type="text" value="8.07"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Reaction Slope	<input type="text" value="+"/>			LIH Influence Check <input type="text" value="No"/>
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	Lipemia <input type="text"/>
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Icterus <input type="text"/>
Linearity Limit	<input type="text"/> %			Hemolysis <input type="text"/>
Lag Time Check	<input type="text"/>			

General	LIH	ISE	Calculated Test	Range				
<b>Test Name:</b> CYS1G <input type="button" value="Test No"/> <b>Type:</b> Serum								
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>					
<b>Specific Ranges</b>								
	Sex	Year	Month	Year	Month	Other Type	Low	High
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7:	Standard demographics						<input type="text" value="#"/>	<input type="text" value="#"/>
8:	Not within expected values						<input type="text" value="#"/>	<input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	<input type="button" value="Select"/>	Decimal Places <input type="text" value="2"/>			

Calibrators	General	ISE				
<b>Test Name:</b> CYS1G <input type="checkbox"/> Use Serum Cal. <b>Type:</b> Serum						
Calibration Type:	<input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>				
<Calibrator Parameters>						
	Calibrator	OD	Conc	Range		Slope Check <input type="text" value="+"/>
				Low	High	
Point-1	CYSC Calibrator Level 1		*	-2.0000	2.0000	Allowable Range Check
Point-2	CYSC Calibrator Level 2		*	-2.0000	2.0000	<input type="checkbox"/> Reagent Blank <input type="text"/>
Point-3	CYSC Calibrator Level 3		*	-2.0000	2.0000	<input type="checkbox"/> Calibration <input type="text"/>
Point-4	CYSC Calibrator Level 4		*	-2.0000	2.0000	Advanced Calibration
Point-5	CYSC Calibrator Level 5		*	-2.0000	2.0000	Operation <input type="text" value="No"/>
Point-6	CYSC Calibrator Level 6		*	-2.0000	2.0000	Interval (RB) <input type="text"/>
Point-7						Interval (ACAL) <input type="text"/>
MB Type Factor	<input type="text"/>	1-Point Calibration Point	<input type="text" value="None"/>	<input type="checkbox"/> with Conc-0		Stability
						Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
						Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

# Definiše korisnik

\* Specifično za seriju, pogledajte tabelu sa analitičkom vrednošću koja je dostupna na [www.gentian.com](http://www.gentian.com)

\*\* Na osnovu rezultata za instrument AU400 (Beckman Coulter®)