Cystatin C

Gentian Cystatin C Immunoassay em sistemas AU da Beckman Coulter® (AU5800, AU680, AU480, DxC 500 AU, DxC 700 AU)

REF B08179

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Referências

1. Abrahamson M et al: Biochem J 1990;268:287-94.
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.
4. Filler G et al: Clin Biochem 2005 ;38 :1-8.
5. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567.
6. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011.
7. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition.* CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute ;2012.
8. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures.* 2nd ed. CLSI guideline EP06. Clinical Laboratory Standards Institut;2020.
9. CLSI. *Establishing and Verifying and Extended Measuring Interval Through Specimen Dilution and Spiking;* 1st Edition. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Institute; 2018.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition.* CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258.
12. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
13. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8.
14. CLSI. *Interference Testing in Clinical Chemistry.* 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
15. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples.* 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
16. Qiu X et al.: Oncotarget. 2017;8(42):72985-72999.
17. CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition.* CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.
18. Nitsch D, et al. Am J Kidney Dis. Jun 2011;57(6):863-72.

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2023-10-12

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www.gentian.com/products/ifu/cystatin-c/beckmancoulter

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Reagente de sistema: B08179 ID de reagente: 228

TEST CONFIGURATION & CHEMISTRY DETAILS																																																	
Assay Name	Test	Rev					Discipline	Chemistry																																									
Test ID	CYS						Calculated Result	<input type="checkbox"/>																																									
LIS Code	CYS																																																
UNITS AND RANGE SETTINGS																																																	
Use Settings from	Serum	Units	mg/L	Decimal Places	x.xx	Plasma																																											
Test Kind	General	Revision	01			<input checked="" type="checkbox"/> Multi Reagent Switch																																											
Reagent Name	CYS	Reagent ID	228			<input type="checkbox"/> FSE Test																																											
ABB Name	CYS1G	Parameter Long Name	Cystatin C B08179 CYS1G CYSC Serum																																														
Region	<input checked="" type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input checked="" type="checkbox"/> AP	<input type="checkbox"/> JP	<input checked="" type="checkbox"/> EU	<input type="checkbox"/> Other																																											
GENERAL PARAMETERS																																																	
SAMPLE VOLUME					REACTION OD LIMIT																																												
Sample Volume	2.0	µL	Dilution	0	µL	Low		High																																									
Predilution Rate	1						REACTION BLANK OD LIMIT																																										
REAGENT VOLUME					ANALYTICAL MEASURING RANGE																																												
R1-1	150	µL	Dilution	0	µL	First: Low	-2.0000	High	2.0000																																								
R2-1	30	µL	Dilution	10	µL	Last: Low	-2.0000	High	2.0000																																								
WAVELENGTH					MANUFACTURER FACTOR																																												
Primary	540	nm	Secondary	None	nm	A	1	B	0																																								
METHOD					REAGENT ONBOARD STABILITY																																												
METHOD					60** Days					0** Hours																																							
REACTION SLOPE					LIH INFLUENCE CHECK																																												
REACTION SLOPE					<input type="checkbox"/> Perform LIH check																																												
MEASURING POINT					Lipemia																																												
Point 1: First	13		Last	27		Icterus	+																																										
Point 2: First			Last			Hemolysis	+																																										
Linearity Limit																																																	
Lag Time Check					<input type="checkbox"/> Perform Lag Time Check																																												
CALIBRATION PARAMETERS																																																	
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4																																								
mg/L	2	None	0	None	0	None	0	None	0																																								
CALIBRATOR SPECIFIC					CALIBRATION OD AND CONCENTRATION PARAMETERS																																												
Calibration Type	6AB	Counts	2	<input type="checkbox"/> Use highest calibrator for Upper AMR																																													
Formula	Spline	MB Factor		<table border="1"> <thead> <tr> <th></th> <th>Calibrator Name</th> <th>Conc</th> <th>OD Range Low</th> <th>OD Range High</th> </tr> </thead> <tbody> <tr> <td>Point 1</td> <td>CYS CAL-1</td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point 2</td> <td>CYS CAL-2</td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point 3</td> <td>CYS CAL-3</td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point 4</td> <td>CYS CAL-4</td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point 5</td> <td>CYS CAL-5</td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point 6</td> <td>CYS CAL-6</td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point 7</td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>							Calibrator Name	Conc	OD Range Low	OD Range High	Point 1	CYS CAL-1	*	-2.0000	2.0000	Point 2	CYS CAL-2	*	-2.0000	2.0000	Point 3	CYS CAL-3	*	-2.0000	2.0000	Point 4	CYS CAL-4	*	-2.0000	2.0000	Point 5	CYS CAL-5	*	-2.0000	2.0000	Point 6	CYS CAL-6	*	-2.0000	2.0000	Point 7				
	Calibrator Name	Conc	OD Range Low	OD Range High																																													
Point 1	CYS CAL-1	*	-2.0000	2.0000																																													
Point 2	CYS CAL-2	*	-2.0000	2.0000																																													
Point 3	CYS CAL-3	*	-2.0000	2.0000																																													
Point 4	CYS CAL-4	*	-2.0000	2.0000																																													
Point 5	CYS CAL-5	*	-2.0000	2.0000																																													
Point 6	CYS CAL-6	*	-2.0000	2.0000																																													
Point 7																																																	
Calibrator Name	CYS	Positive Cutoff																																															
<input checked="" type="checkbox"/> SLOPE CHECK	Slope Check	Number of Levels	6																																														
STABILITY AND INTERVAL					OD DELTA CHECK																																												
Reagent Blank Stability	28 Days	0 Hours	Interval	Lot	<input type="checkbox"/> Reagent Blank	0.0000																																											
Calibration Stability	28 Days	0 Hours	Interval	Lot	<input type="checkbox"/> Calibration	0.0000																																											
PROZONE CHECK PARAMETERS																																																	
<input type="checkbox"/> Logic Check 1			<input type="checkbox"/> Logic Check 2			<input type="checkbox"/> Logic Check 3																																											
Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values																																								
Point 1	#	Value 1	#	Point 1	#	Point 1	#	Point 1	#																																								
Point 2	#	Value 2	#	Interval	#	Value 2	#	Interval	#																																								
Point 3	#	Value 3	#																																														
Limit Points	Limit 1	#	Limit Points	Limit 1	#	Limit Points	Limit 1	Limit 1	#																																								
	Limit 2	#		Limit 2	#		Limit 2	Limit 2	#																																								
Check Pattern	Pattern	#																																															

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Reagente de sistema: B08179 ID de reagente: 228

General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CYS1G"/> <input type="text" value="Test No"/> Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2.0"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit ^{1*} Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> μL	Dilution	<input type="text"/> μL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Analytical Measuring Range Low <input type="text" value="0.40"/> High <input type="text" value="8.07"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/> \downarrow			Onboard Stability Period <input type="text" value="60"/> ** Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/> \downarrow
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CYS1G"/> <input type="text" value="Test No"/> Type: <input type="text" value="Serum"/>				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges				
	Sex <input type="text" value="#"/>	Year <input type="text" value="#"/>	Month <input type="text" value="#"/>	Other Type <input type="text" value="None"/>
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
7:	Standard demographics			<input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	Select <input type="text"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE		
Test Name: <input type="text" value="CYS1G"/> <input type="checkbox"/> Use Serum Cal. Type: <input type="text" value="Serum"/>				
Calibration Type:	<input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>		
<Calibrator Parameters>				
	Calibrator	OD		
		Conc		
		Range		
		Low		
		High		
Point-1	CYSC Calibrator Level 1 \downarrow	*	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2 \downarrow	*	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3 \downarrow	*	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4 \downarrow	*	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5 \downarrow	*	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6 \downarrow	*	-2.0000	2.0000
Point-7	\downarrow			
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/>	<input type="checkbox"/> with Conc-0	Slope Check <input type="text" value="+"/> \downarrow
				Allowable Range Check
				<input type="checkbox"/> Reagent Blank <input type="text"/>
				<input type="checkbox"/> Calibration <input type="text"/>
				Advanced Calibration
				Operation <input type="text" value="No"/> \downarrow
				Interval (RB) <input type="text"/>
				Interval (ACAL) <input type="text"/>
				Stability
				Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
				Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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